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Surface Disinfection Water based

FIELD OF APPLICATION

DCP Surface Disinfection is a rapid and effective water-based disinfectant for use on surfaces. DCP Surface Disinfection is alcohol-free, aldehyde-free, fragrance-free and non-flammable.

INSTRUCTIONS FOR USE

DCP Surface Disinfection is a ready to use, fluid, clear solution. Surfaces and objects are simultaneously und thoroughly moistened, using disposable wipes where appropriate and leaving to dry. Rinsing with water is not necessary. Any remaining product residues are harmless from a toxicological point of view.

SPECIAL APPLICATION

DCP Surface Disinfection can be used in dry fogging applications for surface disinfection.

HIGHLY EFFECTIVE

DCP Surface Disinfection is tested according to European Standards EN 14885, EN 13623, EN 13727, EN 13624, EN 13697, EN 13704, EN 1276 and EN 14476 and is listed by VAH (Verbund für Angewandte Hygiene). DCP Surface Disinfection can be used both as a wipe and spray disinfectant. The stated exposure times below, demonstrate the full effectiveness of DCP Surface Disinfection against bacteria, viruses, spores and yeast.

EXPOSURE TIMES

Bactericidal Effect

precleaned surface

MRSA	1 min	E. coli	1 min*	P. mirabilis.....	1 min*
E. hirae	5 min*	P. aeruginosa	5 min*	L. brevis	5 min
S. aureus	15 min*				

* acc. „Requirements and methods for VAH-Certification (2015/2022)“

Sporocidal Effect

precleaned surface

B. subtilis

5 min

Yeasticidal Effect

precleaned surface

C. albicans

5 min*

T. mentagrophytes ...

5 min

* acc. „Requirements and methods for VAH-Certification (2015/2022)“

Virucidal Effect

Murine Norovirus	clean conditions ... 1 min / dirty conditions ... 1 min
Adenovirus	clean conditions ... 1 min / dirty conditions ... 1 min
Poliovirus	clean conditions ... 1 min / dirty conditions ... 30 min
Vacciniavirus	clean conditions ... 1 min / dirty conditions ... 1 min

Limited Virucidal Effect*

Enveloped viruses incl. human Coronavirus, HIV, HBV, HCV

5 min

(*acc. to the DVV and RKI guidelines)

SPECIFICATIONS

Colour/Form: transparent/liquid pH value: < 7.5

COMPOSITION

>99.5% Water; <0.4% NaCl (Sodium Chloride); <0.055% resp. 0.55g/L HOCl (Hypochlorous acid) / NaOCl (Sodium Hypochlorite)

REGISTRATION BAUA

N-108986 (Product type 2 & 4)

IDENTIFICATION

The product is non-hazardous (classification according (EG) Nr. 1272/2008). Normal precautions must however be observed when handling chemicals. Note the lot number and expiration date printed on the packaging.

SAFETY

In intended use, VDCP Surface Disinfection is not dangerous for people. Avoid eye contact. Irritation may occur. After eye contact, rinse immediately with clear water for several minutes (remove contact lenses). If you experience lengthy complaints related to the product, please consult a medical specialist.

STORAGE

Store DCP Surface Disinfection in a light-protected and tightly closed container. Keep locked up and out of reach of children.

DISPOSAL

With a sufficient dilution DCP Surface Disinfection can be disposed via a standard wastewater sewage system. Please observe local and country specific waste and wastewater regulations. After emptying and thorough rinsing with water, the plastic container can be recycled or can be disposed of in an appropriate way.

SYMBOLS

Temperature limits (3°C – 25°C)

Period after opening in months

Protect from sunlight

Read product information before use

Use disinfectants safely. Always read the label and product information before use.